The background of the slide is a microscopic image of skin tissue, likely a histological section stained with hematoxylin and eosin (H&E). The tissue shows various cellular structures, including what appears to be a basaloid cyst or a similar lesion. Several black circular markers are overlaid on the image, highlighting specific areas of interest. The text is overlaid on this image.

An interoperable, open platform is key for scaled AI implementation in routine workflows

A case study analysis using Proscia DermAI on large datasets

Manuel Garcia Castejon, Director General, Iberia & LATAM | Dr. Thomas Sollie, Medical Director |
Proscia Inc.,
October, 2022

Your speakers for today



Manuel Garcia Castejon

Director General, Iberia & LATAM
Proscia Iberia

Manuel is an industry veteran leading the Spain and LATAM commercial operations for Proscia. Manuel has been a well-known figure in laboratory informatics in the Iberia/ LATAM markets and has helped numerous organizations successfully adopt cutting edge technology solutions to advance lab operations.



Dr. Thomas Sollie

Medical Director

Dr. Thomas Sollie, is the senior clinical AI advisor of pathology for TMC, Unilabs as well as the Medical Director at Proscia Inc., Dr. Thomas Sollie has keen interest in digital and computational pathology and truly believes that technology can help advance cancer diagnostics. Dr. Thomas Sollie brings hands-on experience in implementing digital and computational (AI) pathology solutions at scale across a number of lab sites.



OUR MISSION

Estamos acelerando la transición de la patología a lo digital y utilizando esos datos para cambiar la forma en que diagnosticamos y tratamos el cáncer.

Working globally with our customers to change the way the world practices pathology

Founded in 2014 out of Johns Hopkins



Global reach

Hubs in US and Europe

Nearly 100 Proscians

United in our mission

Experienced leadership

With deep pathology and technology expertise

Domain-focused investors

\$76M+ raised to date

Concentriq digital pathology platform

Frost & Sullivan's 2021 North America Customer Value Leadership Award

Powering the world's digital pathology

Proscia is trusted by **6,500+** pathologists and scientists¹



BIOMARIN®

abbvie



teva



AMGEN



ADVANCED
PATHOLOGY



cernostics
A New View of Cancer Diagnostics



PROTEAN
BIODIAGNOSTICS



UNITY
BIOTECHNOLOGY



UF UNIVERSITY of
FLORIDA



And many more...

¹ 1 June 2021 total registered users across Concentriq Dx and Concentriq for Research

Leading labs trust Proscia's solutions



**First in the world to
100% digital pathology
sign-outs**

“We rigorously evaluated multiple enterprise imaging platforms and chose Concentriq because it gives us full flexibility and infinite scalability with its truly open approach.”



**World's largest repository
of tissue samples
(55 million images)**

JPC is digitizing its tissue archive - the world's largest repository of human pathology specimen - to capitalize on this invaluable source of medical data.



**One of the largest
diagnostic providers in EU
(>200 labs, 200M tests/yr)**

Unilabs has selected Proscia to validate and implement AI-enabled digital pathology (DermAI) in high throughput clinical workflows.

Unified Platform, Purpose-Built Products

Concentriq Digital
Pathology Platform

CONCENTRIQ Dx

Drive efficiency, quality and
confidence in your lab

CONCENTRIQ FOR RESEARCH

Unify research,
accelerate discovery

Concentriq Dx (v 2.1) is CE-marked for in-vitro diagnostic use in Europe and is available for primary diagnosis in the US during the COVID-19 public health emergency.

Construido para potenciar los laboratorios de diagnóstico modernos




Centrado en el patólogo

Una experiencia de usuario centrada en el patólogo




Abierta y verdaderamente interoperable

Interoperabilidad líder con un motor de integración completo



Alta Escalabilidad

Escalabilidad de nivel empresarial con un espacio de trabajo unificado para equipos distribuidos




Computacional de forma nativa

Listo para la inteligencia artificial y entregar patología computacional hoy




Verdaderamente interoperable

Concentriq® Dx se creó para conectar todo su ecosistema de patología digital a través de su API RESTful.




Concentriq Dx (v 2.1) is CE-marked for in-vitro diagnostic use in Europe and is available for primary diagnosis in the US during the COVID-19 public health emergency.




Computacional de forma nativa

La plataforma de patología digital Concentriq® integra la IA directamente en el flujo de trabajo de rutina de alto rendimiento*. La plataforma brinda la oportunidad de ejecutar múltiples aplicaciones de IA automáticamente en función de los metadatos relevantes.

***For research use only**



Concentriq Dx supports scaled interoperability with own & 3P apps




Concentriq Dx (v 4.1.1) is CE-IVDR marked for in-vitro diagnostic use in Europe and is available for primary diagnosis in the US during the COVID-19 public health emergency. Please note that the Galen solution needs to be used in accordance with the intended use and regulatory labelling on the Galen solution as specified by Ibex Medical Analytics



Use of AI in routine practice


Concentriq® Dx and DermAI from Proscia®



Study coverage
 - 134 countries & territories
 - 95% of world population

Data sources
 - national registries
 - publications
 - professional societies
 - personal communication
 - WHO (health workforce)
 - UN (world population)

Global data
 - 13.2M medical doctors
 - 102K+ pathologists
 Asia ≈45K
 N Am ≈25K
 Europe ≈23K
 S Am ≈7K
 Africa <3K
 - 1 pathologist : 125 MD
 US/CAN [1 : 50]
 UK [1 : 100]
 EU [1 : 200]




TOP-10 (70% of workforce)

country	pathologists	per 1M
1. USA	21200+	65
2. India	20000+	14
3. China	11000	8
4. UK	3900	57
5. Brazil	3500	16
6. Russia	3300	73
7. Japan	2600+	21
8. France	2500	29
9. Turkey	2000+	24
10. Australia	1900+	77

data in (no.)
 - definition (ex. Cr in/out)
 - year
 - retired (in/out)
 - residents (in/out)
 - quality of local registries

discussion points
 - lack vs. sufficient
 - workload (volume & pattern)
 - exchange rate
 - income
 - migration

The shortage of pathologists is acute and will become worse!



Technology can be our friend in helping us battle this!

Let us discuss the learnings from the test evaluation of DermAI in a large commercial reference lab setting!



PROSCIA's

Derm AI[®]

Workflow solution designed to drive efficiency in the routine practice of dermatopathology

DermAI is a workflow AI solution that categorizes images of H&E stained Derm specimen into clinically relevant classes - Suspicious Melanocytic, Malignant Basaloid, Malignant Squamous, Benign Melanocytic, Other.

Deployed seamlessly on the CE-IVDR marked Concentriq®Dx image management system. Enables first-read or post-review case sorting, triaging, & prioritization.

It highlights the promise of AI to deliver value in routine practice by:

- Ensuring that cases are routed to the right pathologist
- Helping pathologists prioritize their worklists
- Providing a safety net (AI based second read) on every case for extra confidence

*Proscia DermAI™ is for Research Use Only.

DermAI: One of the most complex AI problems to solve

Dermatopathology is high volume, accounting for **25-30% of case volume** in the lab


Dermatopathology is highly repetitive, with **10 diagnostic outcomes** accounting for **66% of cases**

Melanomas are difficult to read and have **low concordance rates (30-65%)***

*Real world data from DermAI deployment

Distribution of cases is **poorly optimized between generalists and specialists.**

DermAI is built with data from leading labs globally



Trained and validated
on 50,000+ WSIs

Robustly tested for
scanner & stain
variability across global
sites.

PHILIPS

Leica




HAMAMATSU



A case study of DermAI performance on real world data

Test setup (minimal calibration) at one of the largest commercial lab networks

- **Calibration data** set of only 198 WSIs – ground truth confirmed and representative WSIs (1-2) chosen by a Derm pathologist from a total of 195 specimens.
 - An equal number of specimens from every supported DermAI category was used for calibration; namely, 33 melanocytic high-risk specimens, 33 melanocytic low risk specimens, etc.
- **Test data:** A specialist dermatopathologist reviewed an additional 1500 specimens of which WSIs to include and created a **test set with 1695 WSIs**.



For Research Use Only in the United States; requires FDA clearance.

*For research use only



DermAI Performance At-A-Glance

Key benefits observed from DermAI at this lab



30-35% average turnaround time reduction for melanoma (in total)

100% of melanoma cases go directly to Derm paths; 40+% increase in easy cases to general paths



35-40% reduction in number of cases that get a second review.


80% of melanoma cases detected by reviewing only ~16% cases; 70% with just 9% cases.



DermAI is highly accurate in identifying High Risk and Intermediate Risk Melanoma cases

Profile of (GT) Melanoma - High Risk cases as classified by DermAI; All assigned to Dermpaths


Ground Truth Melanocytic High Risk



94.5% of all High-Risk Melanoma cases are classified as High Risk or Suspect by DermAI algorithm

Profile of (GT) Melanoma - Intermediate Risk cases as classified by DermAI; All assigned to Dermpaths

Ground Truth Intermediate Risk



92.7% of all Intermediate-Risk Melanoma cases are classified as High Risk or Suspect by DermAI algorithm




DermAI is also highly specific for other indications as shown below

Category	Sensitivity	Specificity	%Unclassified
Unclassified			5.1%
Melanocytic - High Risk	92%	83%	
Melanocytic - Suspect	91%	83%	
Melanocytic	50%	99%	
BCC	81%	98%	
SCC	76%	92%	




DermAI Performance is proven to be similar across lab and image format variations

Sensitivity, across labs & classification




Specificity, across labs & classification



How useful is DermAI in identifying Melanoma cases?

By looking at 10-15% of incoming cases first, most of the Melanoma cases can be identified









Note: this is a conceptual chart from another US lab site, but we can quickly update this with the actual chart for Unilabs data


- All specimens get a melanoma confidence score from DermAI
- By arranging the cases in the descending order of the Melanoma “confidence score” we identify that:
 - 80% of melanoma cases detected by reviewing only ~16% cases
 - 70% with just 9% cases
- Usually this number is higher (90% melanoma cases identified with ~10% of cases reviewed) in most other laboratories
- This particular data set from this commercial lab was enriched for a higher proportion of Melanoma, so when we test DermAI against a sequential set of cases we expect this performance to be even better as evidences in other laboratories.

Secondly, once these high-risk cases are identified they are re-assigned to specialist dermatopathologists

Without AI
With AI



More importantly, at this lab, DermAI was proven to result in estimated TAT reductions of about 30-35%



>>> 30%+ potential reduction in TAT in most workflows measured till date!


>>> DermAI as a triaging tool is likely ~5X more effective than Clinician flagging

Workflow Impact of DermAI (4)

Estimated Case load distribution & change by class post DermAI for this laboratory


Derm Pathologists
(42% of case load assigned)

Case Load for Dermatopathologists



General Pathologists
(58% of case load assigned)

Case Load for General Pathologists





Key Takeaways

Key Takeaways

1. DermAI calibration burden is minimal and easy to accomplish.
2. DermAI performance is well tested against multiple labs, scanners, and variations. Performance is fairly standard across variations
3. 80% of melanoma cases detected by reviewing only ~16% cases (70% with just 9% of cases reviewed)
4. Estimated 35-40% reduction in cases that are be reviewed twice (first by pathologist and then by dermatopathologist).
5. Estimated 30-35% potential reduction in TAT when DermAI is used in combination with Clinician flagging
6. DermAI as a triaging tool is likely ~5x more effective than Clinician flagging

Thank You



Manuel Garcia Castejon
Director General, Iberia & LATAM, Proscia Inc.,
manuel.co@proscia.com



Dr. Thomas Sollie
Medical Director, Proscia Inc.,
thomas.sollie@proscia.com

Ready for Now. Ready for Next.

proscia.com

